



DRD GUI 6-pack FAVORITES with Dual Res Syringe
DRS 20/300 in PSD/4
Automated PCR FAVORITES
As of 2/11/2008

DRS 20/300 PSD/4 automated PCR Favorites

Edit

Other Favorites:

<p>PCR CP829-ALS</p>  <p>2 5 ul rgnts & 250 nL smpl from 96 plate are aspirated in Diff Mode and delivered together into 1536 plate</p>	<p>PCR CP831-10ALS</p>  <p>2 5 ul rgnts & 250 nL smpl from 96 plate are aspirated Differentially & mixed in tip & delivered to thin 384 PCR</p>	